


*Crystal Systems Corp.*

# Rod Making Package

Model: RM-001



## Crystal Systems Corporation

---

9633-1 Kobuchisawa, Hokuto, Yamanashi 408-0044 JAPAN

TEL +81-551-36-5271 FAX +81-551-36-5273

<http://www.crystalsys.co.jp> [csc@crystalsys.co.jp](mailto:csc@crystalsys.co.jp)

## 1. Outline

This equipment is to help you to make the starting sample rod for Optical Floating Zone method. You can prepare the sample rod utilizing the below-mentioned various parts.

## 2. Parts list

	name	quantity	remarks
1	Vacuum pump (with mouthpiece)	1 unit	
2	Funnel	1 pc	
3	Rubber tubes	100 pcs	$\phi$ 6mm L250mm (one end open)
		10 pcs each	$\phi$ 6,8,10mm L300mm (both ends open)
4	spoon	1 pc	
5	Plastic bar	1 pc	
6	Plastic plate	1 pc	
7	Rubber ring	20pcs	
8	Cotton	1pc	

※ Above content may be changed without prior announcement.